

# Portable Pressure Steam Sterilizer

Electric or LPG heated



## Use

The product utilizes saturated steam to rapidly and efficiently sterilize articles. It can be used in hospitals, public health centers, health stations, clinics in factories and mines, scientific research institutes, etc, to sterilize medical equipment, surgical dressings, glass utensils, solutions, substrata. It can also be used as cooker in plateau areas and be used to make high quality drinking water in industrial and mining enterprises.

## Characteristics

- ◆ Fully stainless steel structure.
- ◆ Electric or LPG heated.
- ◆ Working pressure: 0.14~0.16Mpa.
- ◆ Working temperature: 126℃.
- ◆ Double scale indication pressure gauge.
- ◆ Easy to operate, safe and reliable.

## Technical Data

MODEL	YX-280B18	YX-280B24
Technical Data		
Sterilizing volume	18L φ 280 × 260mm	24L φ 280 × 390mm
Working pressure	0.14~0.16Mpa	
Working temperature	126℃	
Max. Safety pressure	0.165Mpa	
Power	AC220V.50Hz / 2KW	
Transport dimension	410 × 410 × 430	410 × 410 × 550
G.W/N.W	16 / 14kg	17 / 15kg